Public authorities and transfusion medicine are highly concerned by the transfusion Security. KBS created and developed HEMOVIGIL, an intelligent active electronic tag capable of measuring and recording continuously the temperature of labile blood products, from the donor to the final transfusion.

HEMOVIGIL memorizes all the informations related to the blood taken from the donor, its preparation and conservation in the blood bag. This state of the art tracer provides to the practitioner a quality transfusionnal therapy.

By default, HEMOVIGIL monitors the standard parameters required by the legislations concerning the conservation and the transportation of blood bags as prescribes in the Journal Officiel (France), 2004/33/CE and the Guide of good transfusional practices.

HEMOVIGIL continuously indicates the state of conservation of the blood bags thanks to flashing LEDS and without any external device.

The HEMOVIGIL operating system enables the storage of digital informations such as ean-code, transport documents, images, fingerprints, security codes.

The IrDA version of the HEMOVIGIL can be activated by a Pocket PC or a table reader. The collected informations can be introduced into a SQL database through the HTML graphical interface. The traceability informations are generated as pdf files.

The HEMOVIGIL tracer is a real concentrate of technologies, including a high level of technicity, associated with well adapted strategic choices, and providing the most accurate metrological security and stability in the time.

Its internal power of calculation gives him characteristics of analyse, of recording and of transfer unmatched at this day for this kind of product.

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
The CGR ALGORITHM

The life cycle of a blood sample is characterised by these phases : a cooling phase from 38°C to 2-6°C within maximum 24 hours, a storage phase of maximum 41 days (in France) at 2-6°C and a distribution phase where the blood is gradually warmed up from 2-6°C to 24°C within maximum 6 hours

When the proper conditions are fulfilled, the HEMOVIGIL will display a green flashing light, while a red light will be activated when the recommended temperature and time ranges will not be respected

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• Practically, the HEMOVIGIL sticker will be placed on a bag at the moment blood is taken from the donor. Its activation is done by a Pocket PC or an IrDA table interface (all the other starting interventions of the different stages will take place in the same way).

All data concerning the preparation and the processing of the bags are registered in the memory of the HEMOVIGIL

• At the moment the temperature reaches 6°C, the tracer switched automatically in the conservation stage. Inside this conservation stage, one or more transportation stages are possible. The times of departure and arrival will be recorded.

• After the conservation and transportation stages starts the distribution stage, limited to 6 hours and including an automatic programming for the warm-up of the blood to maximum 24 °C.

• Inside the distribution stage is included the effective transfusion. Its time will be recorded together with the numeric informations about the receiver, coming from an internal database.

• After the transfusion, HEMOVIGIL's memory will be downloaded by IrDA and introduced on the web through WIFI or GPRS into an SQL database.

• A stage of reintegration of the bag is included in the HEMOVIGIL allowing, under control and particular validation, to a bag which LED is flashing in red, to be put back again in the circuit. In this case the orange LED will be flashing.

• The biologic data of traceability relative to the contents of the blood bag (taking, laboratory, treatment and qualification) are integrated into HemoVigil in the same way as parameters of calibration of the tracer.

Características técnicas del trazador HemoVigil. Technical characteristics of the HemoVigil.

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